Py4ET 2013: Creating Fixation Density Plot using Heat Maps

Setup Code

• Load necessary Python Modules:

```
In [78]: import numpy as np import matplotlib.cm as cm import matplotlib.pyplot as plt import matplotlib.image as mpimg import brewer2mpl
```

• Switch working directory:

```
In [79]: cd D:\Dropbox\WinPython-32bit-2.7.5.0\my-code\ipython-notebooks
D:\Dropbox\WinPython-32bit-2.7.5.0\my-code\ipython-notebooks
```


Create 2D Gaussian Mask Template

```
In [80]: # Gaussian Mask for Fixations
# Sigma of Gaussian in pixels.
# 33 pixels = ~ 1 visual degree on a 1024x768 monitor at 60 cm.
#
sigma_x = sigma_y = 33.0

# Create x and y pixel ranges for Gauss Mask.
#
x = np.arange(-sigma_x*2.5, sigma_x*2.5, 1)
y = np.arange(-sigma_y*2.5, sigma_y*2.5, 1)

# Create X and Y pixel position values for each element of Gauss. Mask.
#
X, Y = np.meshgrid(x, y)

# Create 2D Gauss Mask as numpy array using X and Y mesh grid data
# and sigma's, with Gauss centered in 2D array (0,0)
# gauss=mlab.bivariate_normal(X, Y, sigma_x, sigma_y, 0,0)
gauss*=1.0/gauss.max()
ghw,ghh=gauss.shape[1]/2,gauss.shape[0]/2
# Lets plot the Gauss Mask for illustrative purposes...
# plt.figure()
gmp=plt.imshow(gauss,cmap=cm.gray,origin='lower',extent=(-ghw, ghw,-ghh, ghh))
```



Load Background Image Displayed During Eye Data Collection

1 of 3 24/07/2013 5:11 AM

```
In [81]: # Load example background image
# Flip vertically...
image_array=numpy.flipud(mpimg.imread("./images/nature_15.png"))

# Get background image size
image_size=(image_array.shape[1],image_array.shape[0])
ihw,ihh=image_size[0]/2,image_size[1]/2

# Display image for illustrative purposes, 0,0 is image center.
#
plt.figure()
bip=plt.imshow(image_array,origin='lower',extent=(-ihw, ihw,-ihh, ihh))
```



Create some Random Fixation Data

```
In [82]: # Here, the fixation event data is being simulated as 500 fixations
                                    # of random position within center 50% of fixation_map,
                                     # and with random fixation durations between 150 and 1500 msec.
                                    sim_fix_count=500
                                     fix_duration_range=150,1500
                                    fixation_x_range=image_size[0]/2, image_size[0]+image_size[0]/2 fixation_y_range=image_size[1]/2, image_size[1]+image_size[1]/2
                                     # Create the dummy Fixation Data as a 3x500 numpy array
                                    fix\_pos=np.column\_stack(\ (np.random.randint(*fixation\_x\_range,\ size=sim\_fix\_count), np.random.randint(*fixation\_y\_range,\ size=sim\_fix\_count), np.random.r
                                    # Plot the fixation data for illustrative purposes.
                                    x = [0,2,4,6,8,10];
                                    y = [0]*len(x);
                                    s = [20*2**n for n in range(len(x))];
                                     plt.figure()
                                     sfd=plt.scatter(fix_pos[:,0],fix_pos[:,1],s=fix_pos[:,2]/10)
                                   1100
                                   1000
                                       900
                                       800
                                       700
                                       600
                                       500
                                       400
```

Create the Fixation Density Map Layer based on the Gauss Mask Template and the Fixation Data

2 of 3 24/07/2013 5:11 AM

```
In [83]: # Start with empty 2D numpy array 2x size of fixation data
    # (makes applying Gauss mask for each fixation easier)
    #
    fixation_map=np.zeros((image_size[1]*2, image_size[0]*2))


# Apply Gauss. Mask to each fixation position based on the created fixation event data.

#
    for fx,fy,fix_dur in fix_pos:
        fixation_map[fy-ghh:fy+ghh+1,fx-ghw:fx+ghh+1]+=(gauss*fix_dur)

# Trim the Fixation Density Mask Back Down to Match the Image Size
    #
    fixation_map=fixation_map[image_size[1]/2:image_size[1]+image_size[1]/2,image_size[0]/2:image_size[0]+image_size[0]/2]


# Plot Fixation Density Mask Alone for illustrative purposes
    # using a Yellow->Orange->Red Color Map.
    clmap=brewer2mpl.get_map('YlOrRd', 'sequential', 8).mpl_colormap
    im = plt.imshow(fixation_map, interpolation='nearest', origin='lower',extent=(-image_size[0]/2, image_size[0]/2,-image_size[1]/2, image_size[1]
    plt.colorbar()
```

Out[83]: <matplotlib.colorbar.Colorbar instance at 0x63639710>



Putting it all Together: Heat Map Representation of Fixation Position and Dwell Time Density During Image Viewing

Out[84]: <matplotlib.colorbar.Colorbar instance at 0x62F357D8>



That's it!

3 of 3 24/07/2013 5:11 AM